

DIAGNOSTIC ABR

IN THE PALM OF YOUR HAND!

**HYBRID MODE!
OPERATE ABR
FROM YOUR PC
OR HANDHELD**



**PATH MEDICAL
celebrates its
10th anniversary!**

**Check out the
offers on our
webpage!**



Connect SENTIERO to your PC via USB and use the new remote control software to display all results live on a large display screen. Use your mouse and QWERTY keyboard to enter patient data and edit protocols in seconds... All diagnostic features are accessible whether connected to a PC or from the handheld device standalone, offering maximum flexibility! Remote control ALL of your test modules, even ASSR! Leave out all the troublesome connectors with just one USB cable required to install SENTIERO into a sound-proof booth. SENTIERO still offers outstanding binaural OAE measurements and the patented DPOAE threshold estimation software, providing an objective cochlear audiogram. Pure tone and speech audiometry tests are extended too! Check out our new speech test modules and live speech option / recorded CD option.

PATH MEDICAL GmbH
Landsberger Straße 65
82110 Germering / Germany

info@pathme.de
www.pathme.de

THE SOUND OF SCIENCE.

